Chemfort™ Closed System Transfer Device (CSTD) Reconstitution and Administration Components

www.simplivia.com

3961555 VER.2020-07

VIAL ADAPTORS (VA) AND CONVERTOR 13mm

DRUG	DILUENT/SALINE			

13mm	20mm	28mm	32mm
VA 20mm converter 13mm	VA 20mm	VA 28mm	VA 32mm

SYRINGE ADAPTOR (SA) **SYRINGE ADAPTOR LOCK (SAL)** **LUER LOCK ADAPTOR (LLA)** **BAG ADAPTOR SP (BASP)** **IV SAFETY SETS**

VIAL ADAPTOR CONNECTION
Each Vial Adaptor fits a specific vial neck

A

2a 20mmVA / 28mmVA / 32mmVA
20mm / 28mm / 32mm

3a CLICK!

2b 13mm
CLICK!

3b 13mm

4b 13mm
CLICK!

SYRINGE ADAPTOR CONNECTION

B

1 MALE LUER LOCK (MLL)

2a

2b

BAG ADAPTOR SP IV SAFETY SETS

C

1

2a

3a

2b

3b Prime

4b Close

DRUG RECONSTITUTION (powder for injection)

D

1a 13mm / 20mm / 28mm / 32mm

2a CLICK!

3a

1b 5 sec IPA 70%

2b

3b

4 CLICK!

5

6

7

DRUG TRANSFER

E

1 13mm / 20mm / 28mm / 32mm

2 CLICK!

3

4

5 AIR ONLY (green checkmark) / No liquid upwards (red X)

6 5 sec IPA 70%

7 5 sec IPA 70%

IV BAG INFUSION

F

1 5 sec IPA 70%

2 CLICK!

3 2-3ml

4

5

6 OPTIONAL: MLL, SA/SAL

ADMINISTRATION-IV DRIP

G

1a **2a**

OR

1b MLL, FEMALE LUER LOCK (FLL)

2b

OR

1c SA/SAL

2c LLA, FLL

3c 5 sec IPA 70%, CLICK!

4c

ADMINISTRATION- IV PUSH/BOLUS

H

1 FLL

2 LLA, FLL

3 CLICK!

4 5 sec IPA 70%

5

GB INDICATIONS FOR USE

Chemfort[™] is a Closed System Transfer Device (CSTD) that mechanically prohibits the release of drugs, including antineoplastic and hazardous drugs, in vapor, aerosol or liquid form during preparation, reconstitution, compounding and administration, minimizing exposure of individuals, healthcare personnel, and the environment to hazardous drugs.

Chemfort[™] prevents the introduction of microbial and airborne contaminants into the drug or fluid path for up to 7 days.

NOTES

- Use Aseptic technique.
- Chemfort**[™] has been shown to prevent microbial ingress for 7 days and 10 activations.
- Replace tubing sets per European Centre for Disease Prevention and Control guidelines/ local regulation and institutional protocol.
- Disposal of components should be according to local regulations and procedures for handling of hazardous drugs.
- Refer to local distributors for additional training information.
- Not made with natural rubber latex or DEHP.

PRECAUTIONS

- Do not disconnect **Chemfort**[™]/**Tevadaptor**[®] devices from other non **Chemfort**[™]/**Tevadaptor**[®] components.
- Refer to drug manufacturer’s and European Pharmacopoeia guidelines for shelf life and sterility information.
- Chemfort**[™] can be used with the drug up to 7 days.

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Simplivia Healthcare Ltd.
Consulting Srl